

MAST Soft X-Ray Cameras

March 27, 2015

Contents

| | | |
|----------|----------------------------------|-----------|
| 1 | Introduction | 2 |
| 1.1 | References | 2 |
| 2 | Configuration | 3 |
| 2.1 | XML Configuration File | 3 |
| 2.2 | C Binding | 3 |
| 2.3 | IDL Binding | 3 |
| 3 | Upper Horizontal Camera | 4 |
| 3.1 | Camera Details | 5 |
| 3.2 | Channel Details | 5 |
| 4 | Lower Horizontal Camera | 6 |
| 4.1 | Camera Details | 6 |
| 4.2 | Channel Details | 7 |
| 5 | Inner Vertical Camera | 8 |
| 5.1 | Camera Details | 8 |
| 5.2 | Channel Details | 9 |
| 6 | Outer Vertical Camera | 9 |
| 6.1 | Camera Details | 9 |
| 6.2 | Channel Details | 10 |
| 7 | Third Horizontal Camera | 10 |
| 7.1 | Camera Details | 11 |
| 7.2 | Channel Details | 11 |
| 8 | Tangential Camera | 12 |
| 8.1 | Camera Details | 12 |
| 8.2 | Channel Details | 13 |

This manual documents the MAST Soft X-Ray Cameras. It is available in [html](http://fuslpcjs.fusion-culham.ukaea.org.uk/sxr/doc/manual/html/index.html) (<http://fuslpcjs.fusion-culham.ukaea.org.uk/sxr/doc/manual/html/index.html>) and pdf formats, generated from xml source. The document tree is stored in a Subversion repository. You can check it out, edit the source and check changes back in, and the html and pdf files will be regenerated. To check out the document tree under Linux do

```
svn co http://fuslpcjs/repos/jstorrs/trunk/sxr/
```

This will create a subdirectory called `sxr` in your current directory, containing the tree. To check in changes to the xml files do

```
svn ci
```

For further information about Subversion see this [tutorial](http://fuslpcjs.fusion.culham.ukaea.org.uk/tutorials/subversion/html/index.html) (<http://fuslpcjs.fusion.culham.ukaea.org.uk/tutorials/subversion/html/index.html>). If you want to generate the html and pdf files in your own workspace, see the [Practical XML](http://fuslpcjs.fusion.culham.ukaea.org.uk/tools/practicalxml/doc/html/index.html) (<http://fuslpcjs.fusion.culham.ukaea.org.uk/tools/practicalxml/doc/html/index.html>) documentation.

Main author: Luca Garzotti.

1 Introduction

Soft X-ray (SXR) diagnostic techniques provide a valuable tool to investigate many properties of plasmas in fusion devices. Applications vary from very simple analysis of line integrated SXR emission allowing for example the determination of the position of the inversion radius associated with saw-tooth crashes to more sophisticated tomographic reconstructions of the shape of the magnetic surfaces inside the plasma.

This document gives details of the SXR cameras on MAST. A schematic of lines of sight and a table of main parameters is provided for each camera. Camera parameters are also available in a C include file and an IDL data file, generated from xml source. They can be used by anyone to develop software tools for the analysis of SXR measurements. The files can be downloaded from the Configuration page.

MAST is equipped with six SXR cameras. Each camera is based on the principle of the pin-hole camera (for a fun introduction to pin-hole cameras see for example [Wikipedia](http://en.wikipedia.org/wiki/Pinhole_camera) (http://en.wikipedia.org/wiki/Pinhole_camera) and links within). It has no focusing optics and is essentially made of a slit and an array of silicon photo-diodes that generate a current when illuminated by photons. The slit is equipped with a beryllium foil to stop the low-energy ($<1\text{keV}$) part of the electromagnetic spectrum and each photo-diode collect the radiation from a specific line of sight across the plasma.

Four of the MAST cameras are situated at sector 12. Two of them sit at the mid-plane port (upper and lower horizontal cameras), one is at the inner top port (outer vertical camera) and another at the outer top port (inner vertical camera). All these cameras have lines of sight lying on the poloidal cross-section of the plasma.

Two more camera are situated at the mid-plane port sector 2. One is similar to the two horizontal cameras (third horizontal camera) with lines of sight lying on the poloidal cross section of the plasma, whereas the other has lines of sight lying in the equatorial plane of the plasma (tangential camera). Because of the limited number of acquisition channels available, until now data were collected either from the third horizontal or from the tangential camera. It is in the plans to increase the number of acquisition channels and start to routinely collect data from both cameras at the same time.

1.1 References

A very nice and instructive paper about SXR tomography on JET is:

Granetz and Smeulders, X-ray tomography on JET, Nuclear Fusion 28 457.

More general introduction to the physics of SXR emission from plasmas can be found in:

Hutchinson, Principle of plasma diagnostics (second ed.), section 5.3, Cambridge university press, 2002.

Presentations on the MAST SXR cameras can be found here:

Turri (2004)

Garzotti (2008)

2 Configuration

2.1 XML Configuration File

Camera configuration data is defined centrally in an XML file. This is processed to generate documentation and code files for use in analysis programs.

2.2 C Binding

A C include file is generated from the XML configuration file. You can include it in your program to use the up-to-date structures defined in it. A small C test program is also provided.

2.3 IDL Binding

The IDL data file is called `sxr.dat`. If you want to use the information in this file you will have first to download it by clicking on the link. Since Windows will try to save it as an mpg file named `sxr.mpg`, on saving the file you will have to rename it `sxr.dat`. You can then restore it in you IDL session by typing:

```
IDL> restore, 'sxr.dat'
```

You will then end up with six structures: UH, LH, IV, OV, TH, TA containing the parameters of the upper horizontal, lower horizontal, inner vertical, outer vertical, third horizontal and tangential camera respectively. You can inspect the content of these structures by typing:

```
IDL> help
% At $MAIN$
IV          STRUCT    = -> <Anonymous> Array[1]
LH          STRUCT    = -> <Anonymous> Array[1]
OV          STRUCT    = -> <Anonymous> Array[1]
TA          STRUCT    = -> <Anonymous> Array[1]
TH          STRUCT    = -> <Anonymous> Array[1]
UH          STRUCT    = -> <Anonymous> Array[1]
Compiled Procedures:
    $MAIN$
Compiled Functions:
```

To see the details of each structure (e. g. IV corresponding to the inner vertical camera) type:

```
IDL> more, IV
IV          STRUCT    = -> <Anonymous> Array[1]
  ID        STRING    = 'IV'
  NAME      STRING    = 'Inner Vertical Camera'
```


| | | | |
|----------------|--------|---|-----------------------|
| FIGURE | STRING | = | 'iv.jpg' |
| SECTOR | INT | = | 12 |
| PORT | STRING | = | 'TO' |
| VIEW | STRING | = | 'Poloidal' |
| SLIT_WIDTH | FLOAT | = | 1.000000 |
| SLIT_HEIGHT | FLOAT | = | 3.000000 |
| DIODE_DISTANCE | FLOAT | = | 62.00000 |
| FOIL_MATERIAL | STRING | = | 'Be' |
| FOIL_THICKNESS | FLOAT | = | 15.00000 |
| DIODE_ARRAY | STRING | = | 'Centronics L35-5T' |
| DIODE_COUNT | INT | = | 12 |
| DIODE_WIDTH | FLOAT | = | 0.9600000 |
| DIODE_HEIGHT | FLOAT | = | 4.600000 |
| POLE_R | FLOAT | = | 700.000 |
| POLE_Z | FLOAT | = | 0.000000 |
| COMMENT | STRING | = | 'Camera description.' |
| CHANNELS | STRUCT | = | -> CHANNEL Array[12] |
| ID | INT | = | Array[12] |
| NAME | STRING | = | Array[12] |
| R1 | FLOAT | = | Array[12] |
| Z1 | FLOAT | = | Array[12] |
| R2 | FLOAT | = | Array[12] |
| Z2 | FLOAT | = | Array[12] |
| P | FLOAT | = | Array[12] |
| PHI | FLOAT | = | Array[12] |
| COMMENT | STRING | = | Array[12] |

The content of each field should be self explanatory, but if in doubt you can ask the RO. It is worth mentioning that, for each element of the structure array CHANNELS, (R1,Z1) and (R2,Z2) are the coordinates in the R-Z plane of two points identifying the line of sight of that particular channel, P is the impact parameter with respect to the pole identified by the coordinates (POLE_R,POLE_Z) and PHI is the angle between the Z=0 axis and the line passing through the pole and perpendicular to that line of sight. These quantities are useful for the tomographic inversion of the SXR signals (see paper by Granetz and Smeulders below). Note also that if you want the parameters P and PHI with respect to a different pole you will have to recalculate them starting from (R1,Z1) and (R2,Z2).

For the tangential camera the situation is slightly different, since the line of sight lie on the machine midplane, instead of on a poloidal cross section. For the tangential camera we have defined a reference frame R-Z where the origin is the centre of the machine, the R axis is the axis separating sector three and sector four and the Z axis is the axis separating sector twelve and sector one.

3 Upper Horizontal Camera

Looks at the plasma below the equatorial plane from the low field side of the machine.



3.1 Camera Details

| | |
|----------------|----------------------|
| id | UH |
| sector | 12 |
| port | HM |
| view | Poloidal |
| slit width | 1.0 |
| slit height | 3.0 |
| diode distance | 50.0 |
| foil material | Be |
| foil thickness | 12.5 |
| diode array | Centronics L35-5T |
| diode count | 18 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 700.0 |
| pole Z | 0.0 |


3.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|-------------|----------|----------|----------|----------|------|--------------|---|
| XSX_HCAMU#1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Channel viewing slightly above midplane. |
| XSX_HCAMU#2 | 196.0 | -39.5 | 2076.0 | 0.0 | 28.9 | 4.733 | SXR emissivity - Centrally viewing channel. |
| XSX_HCAMU#3 | 196.0 | -115.7 | 2076.0 | 0.0 | 84.5 | 4.774 | |

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|--------------|----------|----------|----------|----------|-------|--------------|--------------------------|
| XSX_HCAMU#4 | 196.0 | -186.0 | 2076.0 | 0.0 | 135.5 | 4.811 | |
| XSX_HCAMU#5 | 196.0 | -258.0 | 2076.0 | 0.0 | 187.1 | 4.849 | |
| XSX_HCAMU#6 | 196.0 | -331.0 | 2076.0 | 0.0 | 239.1 | 4.887 | |
| XSX_HCAMU#7 | 196.0 | -407.2 | 2076.0 | 0.0 | 291.1 | 4.926 | |
| XSX_HCAMU#8 | 196.0 | -484.7 | 2076.0 | 0.0 | 343.5 | 4.965 | |
| XSX_HCAMU#9 | 196.0 | -564.0 | 2076.0 | 0.0 | 395.4 | 5.004 | |
| XSX_HCAMU#10 | 196.0 | -657.2 | 2076.0 | 0.0 | 454.1 | 5.049 | |
| XSX_HCAMU#11 | 196.0 | -741.0 | 2076.0 | 0.0 | 504.6 | 5.088 | |
| XSX_HCAMU#12 | 196.0 | -827.0 | 2076.0 | 0.0 | 554.1 | 5.127 | |
| XSX_HCAMU#13 | 196.0 | -915.5 | 2076.0 | 0.0 | 602.4 | 5.166 | |
| XSX_HCAMU#14 | 196.0 | -1006.2 | 2076.0 | 0.0 | 649.3 | 5.204 | |
| XSX_HCAMU#15 | 202.7 | -1095.5 | 2076.0 | 0.0 | 694.6 | 5.242 | |
| XSX_HCAMU#16 | 243.0 | -1165.3 | 2076.0 | 0.0 | 738.2 | 5.279 | |
| XSX_HCAMU#17 | 280.0 | -1253.0 | 2076.0 | 0.0 | 787.3 | 5.322 | Extreme lower channel. |
| XSX_HCAMU#18 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Blanked passive channel. |

4 Lower Horizontal Camera

Looks at the plasma above the equatorial plane from the low field side of the machine.



4.1 Camera Details

| | |
|--------|----|
| id | LH |
| sector | 12 |

| | |
|----------------|----------------------|
| port | HM |
| view | Poloidal |
| slit width | 1.0 |
| slit height | 3.0 |
| diode distance | 50.0 |
| foil material | Be |
| foil thickness | 12.5 |
| diode array | Centronics L35-5T |
| diode count | 18 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 700.0 |
| pole Z | 0.0 |


4.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|--------------|----------|----------|----------|----------|-------|--------------|---|
| XSX_HCAML#1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Channel viewing slightly below midplane. |
| XSX_HCAML#2 | 196.0 | 39.5 | 2076.0 | 0.0 | 28.9 | 1.550 | SXR emissivity - Centrally viewing channel. |
| XSX_HCAML#3 | 196.0 | 115.7 | 2076.0 | 0.0 | 84.5 | 1.509 | |
| XSX_HCAML#4 | 196.0 | 186.0 | 2076.0 | 0.0 | 135.5 | 1.472 | |
| XSX_HCAML#5 | 196.0 | 258.0 | 2076.0 | 0.0 | 187.1 | 1.434 | |
| XSX_HCAML#6 | 196.0 | 331.0 | 2076.0 | 0.0 | 239.1 | 1.396 | |
| XSX_HCAML#7 | 196.0 | 407.2 | 2076.0 | 0.0 | 291.3 | 1.358 | |
| XSX_HCAML#8 | 196.0 | 484.7 | 2076.0 | 0.0 | 343.5 | 1.318 | |
| XSX_HCAML#9 | 196.0 | 564.0 | 2076.0 | 0.0 | 395.4 | 1.279 | |
| XSX_HCAML#10 | 196.0 | 657.2 | 2076.0 | 0.0 | 454.1 | 1.235 | |
| XSX_HCAML#11 | 196.0 | 741.0 | 2076.0 | 0.0 | 504.6 | 1.195 | |
| XSX_HCAML#12 | 196.0 | 827.0 | 2076.0 | 0.0 | 554.1 | 1.156 | |
| XSX_HCAML#13 | 196.0 | 915.5 | 2076.0 | 0.0 | 602.4 | 1.118 | |
| XSX_HCAML#14 | 196.0 | 1006.2 | 2076.0 | 0.0 | 649.3 | 1.079 | |
| XSX_HCAML#15 | 202.7 | 1095.5 | 2076.0 | 0.0 | 694.6 | 1.042 | |
| XSX_HCAML#16 | 243.0 | 1165.3 | 2076.0 | 0.0 | 738.2 | 1.005 | |

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|--------------|----------|----------|----------|----------|-------|--------------|--------------------------|
| XSX_HCAML#17 | 280.0 | 1253.0 | 2076.0 | 0.0 | 787.3 | 0.962 | Extreme upper channel |
| XSX_HCAML#18 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Blanked passive channel. |

5 Inner Vertical Camera

Looks at the plasma high field side and core from the top of the machine.



5.1 Camera Details


| | |
|----------------|----------------------|
| id | IV |
| sector | 12 |
| port | TO |
| view | Poloidal |
| slit width | 1.0 |
| slit height | 3.0 |
| diode distance | 62.0 |
| foil material | Be |
| foil thickness | 15.0 |
| diode array | Centronics L35-5T |
| diode count | 12 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 700.0 |
| pole Z | 0.0 |

5.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|----------------|----------|----------|----------|----------|-------|--------------|--------------------|
| XSX_V_STE29_1 | 966.1 | 2068.4 | 207.3 | 0.0 | 462.6 | 2.790 | Innermost channel. |
| XSX_V_STE29_2 | 966.1 | 2068.4 | 278.3 | 0.0 | 400.2 | 2.821 | |
| XSX_V_STE29_3 | 966.1 | 2068.4 | 348.4 | 0.0 | 336.9 | 2.851 | |
| XSX_V_STE29_4 | 966.1 | 2068.4 | 417.6 | 0.0 | 273.0 | 2.882 | |
| XSX_V_STE29_5 | 966.1 | 2068.4 | 485.8 | 0.0 | 208.6 | 2.913 | |
| XSX_V_STE29_6 | 966.1 | 2068.4 | 553.2 | 0.0 | 144.0 | 2.945 | |
| XSX_V_STE29_7 | 966.1 | 2068.4 | 620.4 | 0.0 | 78.5 | 2.976 | |
| XSX_V_STE29_8 | 966.1 | 2068.4 | 686.0 | 0.0 | 13.9 | 3.007 | |
| XSX_V_STE29_9 | 966.1 | 2068.4 | 750.8 | 0.0 | 50.5 | 6.179 | |
| XSX_V_STE29_10 | 966.1 | 2068.4 | 814.8 | 0.0 | 114.5 | 6.210 | |
| XSX_V_STE29_11 | 966.1 | 2068.4 | 877.9 | 0.0 | 177.7 | 6.241 | |
| XSX_V_STE29_12 | 966.1 | 2068.4 | 940.3 | 0.0 | 240.3 | 6.271 | |

6 Outer Vertical Camera

Looks at the plasma low field side edge from the top of the machine.



6.1 Camera Details

| | |
|----------------|----------|
| id | VO |
| sector | 12 |
| port | TC |
| view | Poloidal |
| slit width | 1.0 |
| slit height | 8.0 |
| diode distance | 120.0 |


| | |
|----------------|----------------------|
| foil material | Be |
| foil thickness | 15.0 |
| diode array | Centronics L35-5T |
| diode count | 12 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 700.0 |
| pole Z | 0.0 |

6.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|----------------|----------|----------|----------|----------|-------|--------------|-------------------------|
| XSX_V_STE36_1 | 712.9 | 2071.8 | 1381.3 | 0.0 | 648.4 | 0.312 | Outermost channel. |
| XSX_V_STE36_2 | 712.9 | 2071.8 | 1344.7 | 0.0 | 616.7 | 0.296 | |
| XSX_V_STE36_3 | 712.9 | 2071.8 | 1308.3 | 0.0 | 584.6 | 0.280 | |
| XSX_V_STE36_4 | 712.9 | 2071.8 | 1272.2 | 0.0 | 552.4 | 0.264 | |
| XSX_V_STE36_5 | 712.9 | 2071.8 | 1236.5 | 0.0 | 520.1 | 0.248 | |
| XSX_V_STE36_6 | 712.9 | 2071.8 | 1201.0 | 0.0 | 487.7 | 0.231 | |
| XSX_V_STE36_7 | 712.9 | 2071.8 | 1165.5 | 0.0 | 454.8 | 0.215 | |
| XSX_V_STE36_8 | 712.9 | 2071.8 | 1131.8 | 0.0 | 423.2 | 0.200 | |
| XSX_V_STE36_9 | 712.9 | 2071.8 | 1095.9 | 0.0 | 389.3 | 0.183 | |
| XSX_V_STE36_10 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Blinded by divertor. |
| XSX_V_STE36_11 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Blinded by divertor. |
| XSX_V_STE36_12 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | Blinded by divertor. |

7 Third Horizontal Camera

Looks at the whole plasma from the low field side of the machine. Located at the same toroidal location as the tangential camera.



7.1 Camera Details

| | |
|----------------|----------------------|
| id | HT |
| sector | 2 |
| port | HM |
| view | Poloidal |
| slit width | 1.0 |
| slit height | 3.0 |
| diode distance | 50.0 |
| foil material | Be |
| foil thickness | 15.0 |
| diode array | Centronics L35-5T |
| diode count | 18 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 700.0 |
| pole Z | 0.0 |


7.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|------------|----------|----------|----------|----------|-------|--------------|------------------------|
| XSX_HPZR_1 | 277.0 | -1249.0 | 2264.0 | 0.0 | 832.3 | 5.274 | Extreme lower channel. |
| XSX_HPZR_2 | 222.0 | -1134.0 | 2264.0 | 0.0 | 759.3 | 5.219 | |
| XSX_HPZR_3 | 198.0 | -998.0 | 2264.0 | 0.0 | 680.3 | 5.162 | |
| XSX_HPZR_4 | 198.0 | -846.0 | 2264.0 | 0.0 | 592.7 | 5.101 | |
| XSX_HPZR_5 | 198.0 | -694.0 | 2264.0 | 0.0 | 498.0 | 5.036 | |

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|-------------|----------|----------|----------|----------|-------|--------------|------------------------|
| XSX_HPZR_6 | 198.0 | -542.0 | 2264.0 | 0.0 | 396.9 | 4.969 | |
| XSX_HPZR_7 | 198.0 | -390.0 | 2264.0 | 0.0 | 290.1 | 4.899 | |
| XSX_HPZR_8 | 198.0 | -239.0 | 2264.0 | 0.0 | 179.7 | 4.828 | |
| XSX_HPZR_9 | 198.0 | -87.0 | 2264.0 | 0.0 | 65.8 | 4.754 | |
| XSX_HPZR_10 | 198.0 | 87.0 | 2264.0 | 0.0 | 65.8 | 1.529 | |
| XSX_HPZR_11 | 198.0 | 239.0 | 2264.0 | 0.0 | 179.7 | 1.456 | |
| XSX_HPZR_12 | 198.0 | 390.0 | 2264.0 | 0.0 | 290.1 | 1.384 | |
| XSX_HPZR_13 | 198.0 | 542.0 | 2264.0 | 0.0 | 396.9 | 1.314 | |
| XSX_HPZR_14 | 198.0 | 694.0 | 2264.0 | 0.0 | 498.0 | 1.247 | |
| XSX_HPZR_15 | 198.0 | 846.0 | 2264.0 | 0.0 | 592.7 | 1.182 | |
| XSX_HPZR_16 | 198.0 | 998.0 | 2264.0 | 0.0 | 680.3 | 1.121 | |
| XSX_HPZR_17 | 222.0 | 1134.0 | 2264.0 | 0.0 | 759.3 | 1.064 | |
| XSX_HPZR_18 | 277.0 | 1249.0 | 2264.0 | 0.0 | 832.3 | 1.010 | Extreme upper channel. |

8 Tangential Camera

Looks at the whole plasma tangentially (as interferometer and HOMER camera).



8.1 Camera Details

| | |
|----------------|------------|
| id | TA |
| sector | 2 |
| port | HM |
| view | Tangential |
| slit width | 1.0 |
| slit height | 3.0 |
| diode distance | 55.0 |

| | |
|----------------|----------------------|
| foil material | Be |
| foil thickness | 15.0 |
| diode array | Centronics L35-5T |
| diode count | 18 |
| diode width | 0.96 |
| diode height | 4.60 |
| pole R | 0.0 |
| pole Z | 0.0 |

8.2 Channel Details

| name | r1 mm | z1 mm | r2 mm | z2 mm | p mm | theta rad | comment |
|-------------|----------|----------|----------|----------|--------|--------------|---|
| XSX_TCAM#1 | 1763.0 | 1441.0 | 0.0 | 0.0 | 0.0 | 0.0 | Extreme inner channel (central column). |
| XSX_TCAM#2 | 1763.0 | 1441.0 | 0.0 | 86.49 | 68.6 | 2.226 | |
| XSX_TCAM#3 | 1763.0 | 1441.0 | 0.0 | 171.60 | 139.3 | 2.195 | |
| XSX_TCAM#4 | 1763.0 | 1441.0 | 0.0 | 256.04 | 212.5 | 2.163 | |
| XSX_TCAM#5 | 1763.0 | 1441.0 | 0.0 | 339.81 | 288.2 | 2.129 | |
| XSX_TCAM#6 | 1763.0 | 1441.0 | 0.0 | 422.93 | 366.3 | 2.094 | |
| XSX_TCAM#7 | 1763.0 | 1441.0 | 0.0 | 505.41 | 446.4 | 2.059 | |
| XSX_TCAM#8 | 1763.0 | 1441.0 | 0.0 | 587.24 | 528.5 | 2.022 | |
| XSX_TCAM#9 | 1763.0 | 1441.0 | 0.0 | 668.46 | 612.3 | 1.984 | |
| XSX_TCAM#10 | 1763.0 | 1441.0 | 0.0 | 760.39 | 709.4 | 1.940 | |
| XSX_TCAM#11 | 1763.0 | 1441.0 | 0.0 | 840.27 | 795.4 | 1.899 | |
| XSX_TCAM#12 | 1763.0 | 1441.0 | 0.0 | 919.54 | 881.8 | 1.858 | |
| XSX_TCAM#13 | 1763.0 | 1441.0 | 0.0 | 998.21 | 968.1 | 1.817 | |
| XSX_TCAM#14 | 1763.0 | 1441.0 | 0.0 | 1076.28 | 1054.0 | 1.775 | |
| XSX_TCAM#15 | 1763.0 | 1441.0 | 0.0 | 1153.77 | 1138.8 | 1.732 | |
| XSX_TCAM#16 | 1763.0 | 1441.0 | 0.0 | 1230.68 | 1222.0 | 1.690 | |
| XSX_TCAM#17 | 1763.0 | 1441.0 | 0.0 | 1308.94 | 1305.3 | 1.646 | |
| XSX_TCAM#18 | 1763.0 | 1441.0 | 0.0 | 1382.77 | 1382.0 | 1.604 | Extreme outer channel (plasma edge). |